

RELEASE NOTE

# Docman 10x

Customer Release

Version 2.229.0

**Joanna James, Principal Product Owner**

01/08/2024



# Disclaimer

Copyright © OneAdvanced 2024.

This document contains confidential and / or proprietary information. The content must not be disclosed to third parties without the prior written approval of Advanced Computer Software Group Limited or one of its subsidiaries as appropriate (each referred to as “OneAdvanced”). External recipients may only use the information contained in this document for the purposes of evaluation of the information and entering into discussions with OneAdvanced and for no other purpose.

Whilst OneAdvanced endeavours to ensure that the information in this document is correct and has been prepared in good faith, the information is subject to change and no representation or warranty is given as to the accuracy or completeness of the information. OneAdvanced does not accept any responsibility or liability for errors or omissions or any liability arising out of its use by external recipients or other third parties.

No information set out or referred to in this document shall form the basis of any contract with an external recipient. Any external recipient requiring the provision of software and/or services shall be required to enter into an agreement with OneAdvanced detailing the terms applicable to the supply of such software and/or services and acknowledging that it has not relied on or been induced to enter into such an agreement by any representation or warranty, save as expressly set out in such agreement.

The software (if any) described in this document is supplied under licence and may be used or copied only in accordance with the terms of such a licence. Issue of this document does not entitle an external recipient to access or use the software described or to be granted such a licence.

The development of OneAdvanced software is continuous and the published information may not reflect the current status. Any particular release of the software may not contain all of the facilities described in this document and / or may contain facilities not described in this document.

Advanced Computer Software Group Limited is a company registered in England and Wales with registration number 05965280 whose registered office is at The Mailbox Level 3, 101 Wharfside Street, Birmingham, B1 1RF.

A full list of its trading subsidiaries is available at [www.oneadvanced.com/privacy-policy/](http://www.oneadvanced.com/privacy-policy/)

## Tech Updates



### Update “Datepicker” package to latest version

Internal reference D10-10897

#### **Description**

As part of Security conformance remediation, the Docman 10 application has been updated to use the latest version of Datepicker. so that we receive relevant security updates and support for the package.



### Remove Cross Site Scripting via Free Text Comments in Tasks

Internal reference D10-11576

#### **Description**

As part of Security conformance remediation, cross site scripting in the Free Text comments of Tasks has been removed, so that the applicable removes vulnerability to an XSS style attack.

## What’s New?



### Patient Facing Services Screen Enhancements

Some UX / UI and text changes have been made to the Patient Facing Services (PFS) functionality. These changes aim to improve the user experience and ensure that users are aware of the patient visibility settings when filing or reviewing a document.



### Filing Screen Text Amendments

Internal reference DCMC-4035

#### **Description**

- The “Do No Share Document” tick box option has been renamed to “Do Not Share Document with Docman Share”
- The “Document Security” field which allows users to select whether to “Allow Access” or “Restrict Access” to the document has been renamed to “Document Access”



### Filing Screen New Tool Tips

Internal reference DCMC-4033 / DCMC-4104

#### **Description**

- A new tooltip icon has been created to appear next to the newly renamed “Do Not Share Document with Docman Share” field. On hovering over the tooltip icon, the text reads “Docman Share makes documents available for viewing by other healthcare organisations outside of the patients registered practice”
- A new tooltip icon has been created to appear next to the “Internal Document Access” field. On hovering over the tooltip icon, the text reads “Use this feature to set who can access the document within the practice”



### Filing Screen Additional Messages

Internal reference DCMC-4040

#### **Description**

- A new warning message in red text has been created to appear beneath the “Create Review” dropdown menu. The warning message states “Online visibility to be set within tasks”
- The warning message will not appear if the user selects “No Review Required” from the “Create Review” dropdown menu.
- The warning message will appear when any other option from the “Create Review” dropdown menu is selected



### Tasks Screen Amendments

Internal references DCMC-4041 / DCMC-4104

#### **Description**


- The “Online Services” section of the Review Task has been moved higher up the page so that it appears underneath the “View in EMIS Web” section
- The buttons in the “Online Services” section have been renamed from “Allow viewing of this document” and “Prevent viewing of this document” to “Visible” and “Non-Visible”
- The UI of the buttons in the “Online Services” sections has been updated so that, when a user selects either “Visible” or “Non-Visible” the button they select will turn blue and the alternate option will be greyed out

**Scotland – Cloud Enabler Tool and Orphan Docs Tool Links Added to Docman Apps page**

Internal reference D10-12018

**Description**

Cloud Enabler application and Orphan Document Tool for Scotland NSS use have been added as links on the [Docman Apps](#) page, as per the screen shot below.




**Workflow Template Options Wording**

Internal reference D10-12207

**Description**

On the Workflow List Screens for Workflow templates, the status drop down text "Inactive/Obsolete" is amended to "Deleted/Obsolete" so that the purpose is clearer for users. (The highlighted item would have the text change)





### Date picker defaults to 19nn

Case 23558669 | Internal reference D10-11629

#### **Description**

When using date picker, the default has been changed to 20nn, unless the year input is less than the current year (in which case, 19nn will be used). This will make date selection quicker and more relevant.

## Defect Fixes



### Quick Steps involving multiple Clinical Codes is not displaying correctly

Case 23395217 | Internal reference D10-10146

#### **Description**

When using quickstep with multiple clinical codes on an Active Coding task, the same clinical code is displayed throughout the quickstep, instead of displaying the next sequential code. Please note, this was a visual issue only, the correct clinical codes were added into audit.



### Improve Error Logging on Patient Synchronisation

Case 23524980 | Internal reference D10-11626

#### **Description**

Site reported issue when synchronising EMIS patients into Docman 10, with error "Value cannot be Null". Additional logging has been provided to establish error when synchronising.



### Enforce Clinical Offset Value for Merge Tool

Case 23691755 | Internal reference D10-12188

#### **Description**

A site encountered issues when merging with the incorrect offset number. The Merge tool has been updated so that the Clinical Offset Value is a required value and matches the value in the configuration.

Kind Regards, The Docman GP Team

**docman**  
delivering paper free care



## Powering the world of work

Our business software is the trusted choice for critical sectors, including healthcare, legal services, and education. We keep the world of work moving.

Speak to our expert consultants for personalised advice & recommendations, & get support on the products you are interested in.

Contact us

 +44(0) 330 343 4000

 [www.oneadvanced.com](http://www.oneadvanced.com)

 [hello@oneadvanced.com](mailto:hello@oneadvanced.com)